

订购热线: 4008-898-798

Anti-AP2B1 antibody

Cat. No. ml125533

Package 25 μl/100 μl/200 μl

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-AP2B1 rabbit polyclonal antibody

Applications ELISA, WB

Immunogen Fusion protein of human AP2B1

Reactivity Human, Mouse, Rat

Content 0.48 mg/ml Host species Rabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol AP2B1

Full name adaptor related protein complex 2 beta 1 subunit

Synonyms ADTB2; AP105B; CLAPB1; AP2-BETA

Swissprot P63010

Target Background

The protein encoded by this gene is one of two large chain components of the assembly protein complex 2, which serves to link clathrin to receptors in coated vesicles. The encoded protein is found on the cytoplasmic face of coated vesicles in the plasma membrane. Two transcript variants encoding different isoforms have been found for this gene.



订购热线: 4008-898-798

Applications

Western blotting

Predicted band size:105 kDa

Positive control:Mouse brain tissue lysate

Recommended dilution: 200-1000

Gel: 6%SDS-PAGE Lysate: 40 µg


Lane: Mouse brain tissue lysate

Primary antibody: ml125533(AP2B1 Antibody) at dilution

1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 10 seconds



ELISA

Recommended dilution:5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn