

Anti-HPS3 antibody

Cat. No.	ml125545
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-HPS3 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human HPS3
Reactivity	Human
Content	0.72 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	HPS3
Full name	HPS3, biogenesis of lysosomal organelles complex 2 subunit 1
Synonyms	SUTAL; BLOC2S1
Swissprot	Q969F9

Target Background

This gene encodes a protein containing a potential clathrin-binding motif, consensus dileucine signals, and tyrosine-based sorting signals for targeting to vesicles of lysosomal lineage. The encoded protein may play a role in organelle biogenesis associated with melanosomes, platelet dense granules, and lysosomes. Mutations in this gene are associated with Hermansky-Pudlak syndrome type 3.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 114 kDa

Positive control: HepG2 cell lysate

Recommended dilution: 200-1000

Gel: 6%SDS-PAGE


Lysate: 40 µg

Lane: HepG2 cell lysate

Primary antibody: ml125545(HPS3 Antibody) at dilution 1/250

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 90 seconds



ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn