

# Anti-TPD52L1 antibody

 Cat. No.
 ml162398

 Package
 25 μl/100 μl/200 μl

 Storage
 -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

| Product overview   |   |
|--------------------|---|
| Description        | Anti-TPD52L1 rabbit polyclonal antibody |
| Applications       | ELISA, WB                               |
| Immunogen          | Synthetic peptide of human TPD52L1      |
| Reactivity         | Human                                   |
| Content            | 0.9 mg/ml                               |
| Host species       | Rabbit                                  |
| lg class           | Immunogen-specific rabbit IgG           |
| Purification       | Antigen affinity purification           |
| Target information |   |
| Currence al        |   |

| Target information |  |
|--------------------|--|
| Symbol             |  |
| Full name          |  |
| Synonyms           |  |
| Swissprot          |  |
|                    |  |

TPD52L1 tumor protein D52-like 1 D53; hD53 Q16890

#### **Target Background**


This gene encodes a member of the tumor protein D52 (TPD52) family. The encoded protein contains a coiled-coil domain and may form homo- or hetero-dimer with TPD52 family members. The protein is reported to be involved in cell proliferation and calcium signaling. It also interacts with the mitogen-activated protein kinase kinase kinase 5 (MAP3K5/ASK1) and positively regulates MAP3K5-induced apoptosis.



## 订购热线: 4008-898-798

Applications Western blotting Predicted band size:22 kDa Positive control:293T, 231, PC3, A549 and K562 cells Recommended dilution: 500-2000

Gel: 12%SDS-PAGE Lysate: 40 µg Lane 1-5: 293T cells, 231 cells, PC3 cells, A549 cells, K562 cells Primary antibody: ml162398(TPD52L1 Antibody) at dilution 1/300 Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution Exposure time: 10 seconds



#### ELISA

Recommended dilution:2000-5000

Wit produce 联系电话: 4008-898-798, 021-61725725

### 联系QQ: 2881505695, 2881505696

- 邮箱: mlbio\_cn@yeah.net
- 网址: www.mlbio.cn