

Anti-IQCD antibody

Cat. No.	ml125770
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-IQCD rabbit polyclonal antibody
Applications	ELISA, WB, IHC
Immunogen	Fusion protein of human IQCD
Reactivity	Human, Mouse
Content	1.02 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	IQCD
Full name	IQ motif containing D
Synonyms	DRC10; CFAP84; 4933433C09Rik
Swissprot	Q96DY2

Target Background

IQCD contains 1 IQ domain. The function of the IQCD protein is not known.

订购热线: 4008-898-798


Applications

Immunohistochemistry

Predicted cell location: Cytoplasm

Positive control: Human tonsil

Recommended dilution: 50-300



The image on the left is immunohistochemistry of paraffin-embedded Human tonsil tissue using ml125770(IQCD Antibody) at dilution 1/60, on the right is treated with fusion protein. (Original magnification: $\times 200$)

Western blotting

Predicted band size: 52 kDa

Positive control: TM4 and 293T cell lysates

Recommended dilution: 1000-5000

Gel: 8%SDS-PAGE


Lysate: 40 μ g

Lane 1-2: TM4 and 293T cell lysates

Primary antibody: ml125770(IQCD Antibody) at dilution 1/1000

Secondary antibody: Goat anti rabbit IgG at 1/5000 dilution

Exposure time: 5 seconds



ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn