

订购热线: 4008-898-798

# Anti-MECR antibody

**Cat. No.** ml125816

**Package** 25 μl/100 μl/200 μl

**Storage** -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

**Product overview** 

**Description** Anti-MECR rabbit polyclonal antibody

Applications ELISA, WB, IHC

Immunogen Fusion protein of human MECR

Reactivity Human, Mouse, Rat

Content1.2 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

**Target information** 

Symbol MECR

Full namemitochondrial trans-2-enoyl-CoA reductaseSynonymsETR1; NRBF1; CGI-63; FASN2B; DYTOABG

Swissprot Q9BV79

#### **Target Background**

The protein encoded by this gene is an oxidoreductase that catalyzes the last step in mitochondrial fatty acid synthesis. Defects in this gene are a cause of childhood-onset dystonia and optic atrophy.





订购热线: 4008-898-798

## **Applications**

# **Immunohistochemistry**

Predicted cell location: Nucleus
Positive control: Human breast cancer
Recommended dilution: 50-300





The image on the left is immunohistochemistry of paraffin-embedded Human breast cancer tissue using ml125816(MECR Antibody) at dilution 1/70, on the right is treated with fusion protein. (Original magnification: ×200)

## Western blotting

Predicted band size:40 kDa

Positive control:Human fetal brain tissue lysate

Recommended dilution: 1000-5000

Gel: 8%SDS-PAGE

Lysate: 40 µg

Lane: Human fetal brain tissue lysate

Primary antibody: ml125816(MECR Antibody) at dilution 1/1000

Secondary antibody: Goat anti rabbit IgG at 1/5000 dilution

Exposure time: 10 seconds

kDa
95 —
72 —
55 —
36 —
28 —

### **ELISA**

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio cn@yeah.net

网址: www.mlbio.cn