

订购热线: 4008-898-798

Anti-SMPD2 antibody

Cat. No. ml125865

Package 25 μl/100 μl/200 μl

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-SMPD2 rabbit polyclonal antibody

Applications ELISA, WB, IHC

Immunogen Fusion protein of human SMPD2

Reactivity Human, Mouse, Rat

Content 1.26 mg/ml Host species Rabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol SMPD2

Full name sphingomyelin phosphodiesterase 2

Synonyms ISC1; NSMASE; NSMASE1

Swissprot O60906

Target Background

This gene encodes a protein which was initially identified as a sphingomyelinase based on sequence similarity between bacterial sphingomyelinases and a yeast protein. Subsequent studies showed that its biological function is less likely to be as a sphingomyelinase and instead as a lysophospholipase.





订购热线: 4008-898-798

Applications


Immunohistochemistry

Predicted cell location: Cytoplasm Positive control: Human liver cancer Recommended dilution: 50-300





The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using ml125865(SMPD2 Antibody) at dilution 1/70, on the right is treated with fusion protein. (Original magnification: ×200)



Predicted band size:48 kDa

Positive control: Mouse lung tissue lysate Recommended dilution: 1000-5000

Gel: 8%SDS-PAGE

Lysate: 40 µg


Lane: Mouse lung tissue lysate


on di Primary antibody: ml125865(SMPD2 Antibody) at dilution 1/1000

Secondary antibody: Goat anti rabbit IgG at 1/5000 dilution


Exposure time: 2 minutes

Predicted cell location: Cytoplasm Positive control: Human tonsil Recommended dilution: 50-300





The image on the left is immunohistochemistry of paraffin-embedded Human tonsil tissue using ml125865(SMPD2 Antibody) at dilution 1/70, on the right is treated with fusion protein. (Original magnification: ×200)



ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio cn@yeah.net

网址: www.mlbio.cn