

Anti-XKR3 antibody

Cat. No.	ml162622
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-XKR3 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human XKR3
Reactivity	Human
Content	0.2 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	XKR3
Full name	XK, Kell blood group complex subunit-related family, member 3
Synonyms	XRG3; XTES
Swissprot	Q5GH77

Target Background

XKR3 (MIM 300684) and XKR3 are homologs of the Kell blood group precursor XK (MIM 314850), which is a putative membrane transporter and a component of the XK/Kell complex of the Kell blood group system (Calenda et al., 2006 [PubMed 16431037]).

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 53 kDa

Positive control: Human normal stomach and thyroid cancer tissue

Recommended dilution: 500-2000

Gel: 8%SDS-PAGE


Lysate: 40 µg

Lane 1-2: Human normal stomach and thyroid cancer tissue

Primary antibody: ml162622(XKR3 Antibody) at dilution 1/330

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 2 minutes



ELISA

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn