

Anti-ZFAND2B antibody

Cat. No.	ml162630
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-ZFAND2B rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human ZFAND2B
Reactivity	Human, Mouse, Rat
Content	0.4 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	ZFAND2B
Full name	zinc finger, AN1-type domain 2B
Synonyms	AIRAPL
Swissprot	Q8WV99

Target Background

This gene encodes a protein containing AN1-type zinc-fingers and ubiquitin-interacting motifs. The encoded protein likely associates with the proteasome to stimulate the degradation of toxic or misfolded proteins. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 28 kDa

Positive control: Mouse brain tissue

Recommended dilution: 200-1000

Gel: 12% SDS-PAGE


Lysate: 40 µg

Lane: Mouse brain tissue

Primary antibody: ml162630 (ZFAND2B Antibody) at dilution 1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute



ELISA

Recommended dilution: 1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn