

Anti-GJB5 antibody

Cat. No.	ml162692
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-GJB5 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human GJB5
Reactivity	Human, Mouse, Rat
Content	0.7 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	GJB5
Full name	gap junction protein, beta 5, 31.1kDa
Synonyms	CX31.1
Swissprot	O95377

Target Background

This gene encodes a member of the beta-type (group I) connexin family. The encoded protein is a gap junction protein involved in intercellular communication related to epidermal differentiation and environmental sensing. This gene has been linked to non-small cell lung cancer.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size:31 kDa

Positive control:Human placenta tissue

Recommended dilution: 200-1000

Gel: 10%SDS-PAGE


Lysate: 40 μ g

Lane: Human placenta tissue

Primary antibody: ml162692(GJB5 Antibody) at dilution 1/400

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 15 seconds



ELISA

Recommended dilution:1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn