

Anti-ZNRF1 antibody

Cat. No.	ml152125
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-ZNRF1 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human ZNRF1
Reactivity	Human, Mouse
Content	2.9 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	ZNRF1
Full name	zinc and ring finger 1, E3 ubiquitin protein ligase
Synonyms	NIN283
Swissprot	Q8ND25

Target Background

This gene encodes an E3 ubiquitin-protein ligase that plays a role in neural-cell differentiation. Overexpression of this gene causes neurite-like elongation. The encoded protein contains both a zinc finger and a RING finger motif and is localized in the endosome/lysosome compartment, indicating that it may be involved in ubiquitin-mediated protein modification, and in synaptic vesicle membranes in neurons.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 24 kDa

Positive control: Hela cells

Recommended dilution: 200-1000

Gel: 10% SDS-PAGE


Lysate: 30 µg

Lane: Hela cells

Primary antibody: ml152125(ZNRF1 Antibody) at dilution 1/500

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 10 seconds



ELISA

Recommended dilution: 1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn