

Anti-SLC30A6 antibody

Cat. No.	ml152127
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-SLC30A6 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human SLC30A6
Reactivity	Human, Mouse
Content	1 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	SLC30A6
Full name	solute carrier family 30 (zinc transporter), member 6
Synonyms	ZNT6; MST103; MSTP103
Swissprot	Q6NXT4

Target Background

Zinc functions as a cofactor for numerous enzymes, nuclear factors, and hormones and as an intra- and intercellular signal ion. Members of the zinc transporter (ZNT)/SLC30 subfamily of the cation diffusion facilitator family, such as SLC30A6, permit cellular efflux of zinc.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 51 kDa

Positive control: Mouse thymus tissue

Recommended dilution: 500-2000

Gel: 10% SDS-PAGE


Lysate: 30 µg

Lane: Mouse thymus tissue

Primary antibody: ml152127 (SLC30A6 Antibody) at dilution 1/223

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 2 minutes



ELISA

Recommended dilution: 1000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn