

订购热线: 4008-898-798

Anti-ZBTB10 antibody

Cat. No. ml156008

Package 25 μl/100 μl/200 μl

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-ZBTB10 rabbit polyclonal antibody

Applications ELISA, WB

Immunogen Synthetic peptide of human ZBTB10

ReactivityHuman, RatContent0.5 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol ZBTB10

Full name zinc finger and BTB domain containing 10

Synonyms RINZF Swissprot Q96DT7

Target Background

RINZF, also known as ZBTB10 (zinc finger and BTB domain containing protein 10), is a 847 amino acid protein that contains one BTB/POZ domain and two C2H2-type zinc fingers. Localized to the nucleus, RINZF is believed to play a role in transcriptional regulation. Specifically, RINZF is capable of binding to the CACC element of the Gastrin promoter. In this regard, RINZF competes with Sp1 for CACC binding and interferes with Sp1 transactivation, thereby regulating Gastrin gene expression. The rat RINZF protein shares 98% homology with the human RINZF protein, suggesting that RINZF is a conserved protein. Due to alternative splicing events, two RINZF isoforms exist. In addition, RINZF may be phosphorylated by ATR or ATM upon DNA damage.



订购热线: 4008-898-798

Applications

Western blotting

Predicted band size:95 kDa Positive control:293T, Jurkat, Raji and hela cells

Recommended dilution: 500-2000

Gel: 8%SDS-PAGE

Lysate: 40 µg


Lane 1-4: 293T cells, Jurkat cells, Raji cells, hela cells

Primary antibody: ml156008(ZBTB10 Antibody) at

dilution 1/1000

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 40 seconds



ELISA

Recommended dilution:2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn