

Anti-ZFP28 antibody

Cat. No.	ml156025
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-ZFP28 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human ZFP28
Reactivity	Human, Mouse
Content	0.9 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	ZFP28
Full name	ZFP28 zinc finger protein
Synonyms	mkr5
Swissprot	Q8NHY6

Target Background

May be involved in transcriptional regulation. May have a role in embryonic development.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 99 kDa

Positive control: 293T cells

Recommended dilution: 200-1000

Gel: 8%SDS-PAGE


Lysate: 40 µg

Lane: 293T cells

Primary antibody: ml156025(ZFP28 Antibody) at dilution 1/500

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 2 minutes



ELISA

Recommended dilution: 1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn