

订购热线: 4008-898-798

# Anti-SERPINH1 antibody

**Cat. No.** ml162935

**Package** 25 μl/100 μl/200 μl

**Storage** -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

**Product overview** 

**Description** Anti-SERPINH1 rabbit polyclonal antibody

Applications ELISA, WB

Immunogen Synthetic peptide of human SERPINH1

Reactivity Human, Mouse, Rat

Content0.5 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

**Target information** 

Symbol SERPINH1

Full name serpin peptidase inhibitor, clade H (heat shock protein 47), member 1, (collagen binding

protein 1)

**Synonyms** CBP1; CBP2; Ol10; gp46; AsTP3; HSP47; PIG14; PPROM; RA-A47; SERPINH2

Swissprot P50454

#### **Target Background**

This gene encodes a member of the serpin superfamily of serine proteinase inhibitors. The encoded protein is localized to the endoplasmic reticulum and plays a role in collagen biosynthesis as a collagen-specific molecular chaperone. Autoantibodies to the encoded protein have been found in patients with rheumatoid arthritis. Expression of this gene may be a marker for cancer, and nucleotide polymorphisms in this gene may be associated with preterm birth caused by preterm premature rupture of membranes.



订购热线: 4008-898-798

## **Applications**

## Western blotting

Predicted band size:46 kDa Positive control:NIH/3T3 cell Recommended dilution: 500-2000

Gel: 8%SDS-PAGE Lysate: 40 µg Lane: NIH/3T3 cell

Primary antibody: ml162935(SERPINH1 Antibody) at

dilution 1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 5 seconds

kDa 130— 95— 72— 55— 36— 28— 17—

### **ELISA**

Recommended dilution:5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio\_cn@yeah.net

网址: www.mlbio.cn