

Anti-Gzmd antibody

Cat. No.	ml162941
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-Gzmd rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of mouse Gzmd
Reactivity	Mouse
Content	0.3 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	Gzmd
Full name	granzyme D
Synonyms	CCP5; Ctla5; Ctla-5
Swissprot	P11033

Target Background

This enzyme is probably necessary for target cell lysis in cell-mediated immune responses.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 28 kDa

Positive control: Mouse heart tissue

Recommended dilution: 200-1000

Gel: 12% SDS-PAGE


Lysate: 40 µg

Lane: Mouse heart tissue

Primary antibody: ml162941 (Gzmd Antibody) at dilution 1/350

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 10 seconds



ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn