

Anti-UBXN2A antibody

Cat. No.	ml160074
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-UBXN2A rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human UBXN2A
Reactivity	Human, Mouse
Content	0.5 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	UBXN2A
Full name	UBX domain-containing protein 2A
Synonyms	UBXD4
Swissprot	P68543

Target Background

The protein is located in both the ER and cis-Golgi compartments from biochemical and immunofluorescence studies.UBXD4 as one of the cytosolic proteins that interact directly with the α3 and α4 nAChR subunits.Overexpression of UBXD4 in differentiated PC12 cells (dPC12) increased nAChR cell surface expression, especially that of the α3β2 subtype.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 29 kDa

Positive control: Mouse liver tissue

Recommended dilution: 200-1000

Gel: 10% SDS-PAGE


Lysate: 40 µg

Lane: Mouse liver tissue

Primary antibody: ml160074 (UBXN2A Antibody) at dilution 1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 10 seconds



ELISA

Recommended dilution: 1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn