

Anti-CPT2 antibody

Cat. No.	ml160201
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-CPT2 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human CPT2
Reactivity	Human, Mouse, Rat
Content	0.7 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	CPT2
Full name	Carnitine palmitoyltransferase 2
Synonyms	CPT1, IIAE4, CPTASE
Swissprot	P23786

Target Background

The protein encoded by this gene is a nuclear protein which is transported to the mitochondrial inner membrane. Together with carnitine palmitoyltransferase I, the encoded protein oxidizes long-chain fatty acids in the mitochondria. Defects in this gene are associated with mitochondrial long-chain fatty-acid (LCFA) oxidation disorders.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 74 kDa

Positive control: Mouse kidney tissue

Recommended dilution: 500-2000

Gel: 6% SDS-PAGE


Lysate: 40 µg

Lane: Mouse kidney tissue

Primary antibody: ml160201(CPT2 Antibody) at dilution 1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 15 seconds



ELISA

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn