

Anti-FOSL1 antibody

Cat. No.	ml160339
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-FOSL1 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human FOSL1
Reactivity	Human, Rat
Content	0.7 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	FOSL1
Full name	FOS-like antigen 1
Synonyms	FRA; FRA1; fra-1
Swissprot	P15407

Target Background

The Fos gene family consists of 4 members: FOS, FOSB, FOSL1, and FOSL2. These genes encode leucine zipper proteins that can dimerize with proteins of the JUN family, thereby forming the transcription factor complex AP-1. As such, the FOS proteins have been implicated as regulators of cell proliferation, differentiation, and transformation.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 29 kDa

Positive control: 231 cells

Recommended dilution: 500-2000

Gel: 10% SDS-PAGE


Lysate: 40 µg

Lane: 231 cells

Primary antibody: ml160339 (FOSL1 Antibody) at dilution 1/615

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 20 seconds



ELISA

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn