

## Anti-FOXD3 antibody

<b>Cat. No.</b>	ml160351
<b>Package</b>	25 µl/100 µl/200 µl
<b>Storage</b>	-20°C, pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol

### Product overview

<b>Description</b>	Anti-FOXD3 rabbit polyclonal antibody
<b>Applications</b>	ELISA, WB
<b>Immunogen</b>	Synthetic peptide of human FOXD3
<b>Reactivity</b>	Human, Mouse
<b>Content</b>	0.3 mg/ml
<b>Host species</b>	Rabbit
<b>Ig class</b>	Immunogen-specific rabbit IgG
<b>Purification</b>	Antigen affinity purification

### Target information

<b>Symbol</b>	FOXD3
<b>Full name</b>	forkhead box D3
<b>Synonyms</b>	AIS1, HFH2, VAMAS2, Genesis
<b>Swissprot</b>	Q9UJU5

### Target Background

Forkhead box D3 also known as FOXD3 is a forkhead protein that in humans is encoded by the FOXD3 gene. This gene belongs to the forkhead protein family of transcription factors which is characterized by a DNA-binding forkhead domain. FoxD3 functions as a transcriptional repressor and contains the C-terminal engrailed homology-1 motif (eh1), which provides an interactive surface with a transcriptional co-repressor Grg4 (Groucho-related gene-4).

订购热线: 4008-898-798

**Applications**

**Western blotting**

Predicted band size:48 kDa

Positive control:Hela cells

Recommended dilution: 500-2000

Gel: 8%SDS-PAGE


Lysate: 40  $\mu$ g

Lane: Hela cells

Primary antibody: ml160351(FOXD3 Antibody) at dilution 1/450

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute



**ELISA**

Recommended dilution:1000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio\_cn@yeah.net

网址: www.mlbio.cn