

订购热线: 4008-898-798

# Anti-DACH1 antibody

**Cat. No.** ml122197

**Package** 25 μl/100 μl/200 μl

**Storage** -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

**Product overview** 

**Description** Anti-DACH1 rabbit polyclonal antibody

Applications ELISA, IHC

Immunogen Fusion protein of human DACH1

ReactivityHuman, MouseContent0.8 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

**Target information** 

Symbol DACH1

Full name dachshund homolog 1 (Drosophila)

Synonyms DACH
Swissprot Q9UI36

### **Target Background**

This gene encodes a chromatin-associated protein that associates with other DNA-binding transcription factors to regulate gene expression and cell fate determination during development. The protein contains a Ski domain that is highly conserved from Drosophila to human. Expression of this gene is lost in some forms of metastatic cancer, and is correlated with poor prognosis. Multiple transcript variants encoding different isoforms have been found for this gene.




订购热线: 4008-898-798

## **Applications**

## **Immunohistochemistry**

Predicted cell location: Nucleus Positive control: Human tonsil Recommended dilution: 50-200







The image on the left is immunohistochemistry of paraffin-embedded Human tonsil tissue using ml122197(DACH1 Antibody) at dilution 1/50, on the right is treated with fusion protein. (Original magnification: ×200)

Predicted cell location: Nucleus Positive control: Human thyroid cancer

Recommended dilution: 50-200





The image on the left is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using ml122197(DACH1 Antibody) at dilution 1/50, on the right is treated with fusion protein. (Original magnification: ×200)

#### **ELISA**

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio\_cn@yeah.net 网址: www.mlbio.cn