

订购热线: 4008-898-798

Anti-CELF5 antibody

Cat. No. ml163294

Package 25 μl/100 μl/200 μl

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-CELF5 rabbit polyclonal antibody

Applications ELISA, WB

Immunogen Synthetic peptide of human CELF5

ReactivityHumanContent0.4 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol CELF5

Full nameCUGBP, Elav-like family member 5SynonymsCELF-5; BRUNOL5; BRUNOL-5

Swissprot Q8N6W0

Target Background

This gene encodes a member of the the CELF/BRUNOL protein family, which contain two N-terminal RNA recognition motif (RRM) domains, one C-terminal RRM domain, and a divergent segment of 160-230 aa between the second and third RRM domains. Members of this protein family regulate pre-mRNA alternative splicing and may also be involved in mRNA editing and translation. Alternatively spliced transcript variants have been found for this gene.



订购热线: 4008-898-798

Applications

Western blotting

Predicted band size:52 kDa
Positive control:A172 cell lysate
Recommended dilution: 200-1000

Gel: 8%SDS-PAGE

Lysate: 40 µg


Lane: A172 cell lysate

Primary antibody: ml163294(CELF5 Antibody) at dilution

1/300

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute



ELISA

Recommended dilution:5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn