

Anti-RCAN3 antibody

Cat. No.	ml163363
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-RCAN3 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human RCAN3
Reactivity	Human, Mouse
Content	0.4 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	RCAN3
Full name	RCAN family member 3
Synonyms	RCN3; MCIP3; hRCN3; DSCR1L2
Swissprot	Q9JKK0

Target Background

RCAN3 inhibits calcineurin-dependent transcriptional responses by binding to the catalytic domain of calcineurin A. It could play a role during central nervous system development.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size:27 kDa

Positive control:Mouse brain tissue lysate

Recommended dilution: 200-1000

Gel: 12%SDS-PAGE


Lysate: 40 µg

Lane: Mouse brain tissue lysate

Primary antibody: ml163363(RCAN3 Antibody) at dilution 1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute



ELISA

Recommended dilution:5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn