

Anti-HINT3 antibody

Cat. No.	ml160478
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-HINT3 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human HINT3
Reactivity	Human
Content	0.2 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	HINT3
Full name	histidine triad nucleotide binding protein 3
Synonyms	HINT4
Swissprot	Q9NQE9

Target Background

Histidine triad proteins, such as HINT3, are nucleotide hydrolases and transferases that act on the alpha-phosphate of ribonucleotides. Hydrolyzes phosphoramidate and acyl-adenylate substrates.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 20 kDa

Positive control: Human prostate tissue

Recommended dilution: 200-1000

Gel: 12% SDS-PAGE


Lysate: 40 µg

Lane: Human prostate tissue

Primary antibody: ml160478 (HINT3 Antibody) at dilution 1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 2 minutes



ELISA

Recommended dilution: 1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn