

Anti-DCAF13 antibody

Cat. No.	ml163675
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-DCAF13 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human DCAF13
Reactivity	Human, Mouse
Content	0.66 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	DCAF13
Full name	DDB1 and CUL4 associated factor 13
Synonyms	GM83; WDSOF1; HSPC064
Swissprot	Q9NV06

Target Background

DCAF13 localizes to the granular component of the nucleolus, a region associated with late ribosomal RNA-processing events and assembly of ribosomal particles. Possible role in ribosomal RNA processing (By similarity). May function as a substrate receptor for CUL4-DDB1 E3 ubiquitin-protein ligase complex.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 51 kDa

Positive control: Human fetal liver tissue lysate

Recommended dilution: 200-1000

Gel: 8%SDS-PAGE


Lysate: 40 µg

Lane: Human fetal liver tissue lysate

Primary antibody: ml163675(DCAF13 Antibody) at dilution 1/250

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute



ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn