

Anti-OTUD6A antibody

Cat. No.	ml163743
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-OTUD6A rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human OTUD6A
Reactivity	Human, Mouse
Content	1.44 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	OTUD6A
Full name	OTU deubiquitinase 6A
Synonyms	DUBA2; HSHIN6
Swissprot	Q7L8S5

Target Background

Deubiquitinating enzymes (DUBs; see MIM 603478) are proteases that specifically cleave ubiquitin (MIM 191339) linkages, negating the action of ubiquitin ligases. DUBA2 belongs to a DUB subfamily characterized by an ovarian tumor (OTU) domain.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 33 kDa

Positive control: Human breast cancer tissue, TM4, 231 and Jurkat cell lysates

Recommended dilution: 500-2000

Gel: 8%SDS-PAGE


Lysate: 40 µg

Lane 1-4: Human breast cancer tissue, TM4, 231 and Jurkat cell lysates

Primary antibody: ml163743(OTUD6A Antibody) at dilution 1/550

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 20 seconds



ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn