

Anti-PFN1 antibody

Cat. No.	ml160711
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-PFN1 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human PFN1
Reactivity	Human, Mouse, Rat
Content	0.6 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	PFN1
Full name	profilin 1
Synonyms	ALS18
Swissprot	P07737

Target Background

This gene encodes a member of the profilin family of small actin-binding proteins. The encoded protein plays an important role in actin dynamics by regulating actin polymerization in response to extracellular signals. Deletion of this gene is associated with Miller-Dieker syndrome, and the encoded protein may also play a role in Huntington disease. Multiple pseudogenes of this gene are located on chromosome 1.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 15 kDa

Positive control: Human seminoma tissue, hela and A549 cells

Recommended dilution: 500-2000

Gel: 12%SDS-PAGE


Lysate: 10 µg

Lane 1-3: Human seminoma tissue, hela cells, A549 cells

Primary antibody: ml160711(PFN1 Antibody) at dilution 1/1000

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute



ELISA

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn