

订购热线: 4008-898-798

Anti-TNFRSF17 antibody

Cat. No. ml160956

Package 25 μl/100 μl/200 μl

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-TNFRSF17 rabbit polyclonal antibody

Applications ELISA, WB

Immunogen Synthetic peptide of human TNFRSF17

ReactivityHuman, MouseContent0.6 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol TNFRSF17

Full name tumor necrosis factor receptor superfamily, member 17

Synonyms BCM; BCMA; CD269; TNFRSF13A

Swissprot Q02223

Target Background

The protein encoded by this gene is a member of the TNF-receptor superfamily. This receptor is preferentially expressed in mature B lymphocytes, and may be important for B cell development and autoimmune response. This receptor has been shown to specifically bind to the tumor necrosis factor (ligand) superfamily, member 13b (TNFSF13B/TALL-1/BAFF), and to lead to NF-kappaB and MAPK8/JNK activation. This receptor also binds to various TRAF family members, and thus may transduce signals for cell survival and proliferation



订购热线: 4008-898-798

Applications

Western blotting

Predicted band size:20 kDa

Positive control:Mouse skeletal muscle tissue

Recommended dilution: 200-1000

Gel: 10%SDS-PAGE

Lysate: 40 µg


Lane: Mouse skeletal muscle tissue

Primary antibody: ml160956(TNFRSF17 Antibody) at

dilution 1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 15 seconds



ELISA

Recommended dilution:1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn