

Anti-ARHGAP45 antibody

Cat. No.	ml164144
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-ARHGAP45 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human ARHGAP45
Reactivity	Human
Content	1.56 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	ARHGAP45
Full name	Rho GTPase activating protein 45
Synonyms	HA-1; HMHA1; HLA-HA1
Swissprot	Q92619

Target Background

Contains a GTPase activator for the Rho-type GTPases (RhoGAP) domain that would be able to negatively regulate the actin cytoskeleton as well as cell spreading. However, also contains N-terminally a BAR-domain which is able to play an autoinhibitory effect on this RhoGAP activity.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 125 kDa

Positive control: A375 cell lysate

Recommended dilution: 500-2000

Gel: 6%SDS-PAGE


Lysate: 40 μ g

Lane: A375 cell lysate

Primary antibody: ml164144 (ARHGAP45 Antibody) at dilution 1/700

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute



ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn