

Anti-HOOK3 antibody

Cat. No.	ml164206
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-HOOK3 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human HOOK3
Reactivity	Human, Mouse
Content	0.66 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	HOOK3
Full name	hook microtubule tethering protein 3
Synonyms	HK3
Swissprot	Q86VS8

Target Background

Hook proteins are cytosolic coiled-coil proteins that contain conserved N-terminal domains, which attach to microtubules, and more divergent C-terminal domains, which mediate binding to organelles. The Drosophila Hook protein is a component of the endocytic compartment.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 83 kDa

Positive control: Human cerebrum tissue, Human cerebella tissue lysates

Recommended dilution: 200-1000

Gel: 6%SDS-PAGE


Lysate: 40 µg

Lane 1-2: Human cerebrum tissue, Human cerebella tissue lysates

Primary antibody: ml164206(HOOK3 Antibody) at dilution 1/300

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute



ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn