

Anti-IFNA5 antibody

Cat. No.	ml164234
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-IFNA5 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human IFNA5
Reactivity	Human
Content	1.56 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	IFNA5
Full name	interferon alpha 5
Synonyms	INA5; INFA5; IeIF G; IFN-alphaG; IFN-alpha-5
Swissprot	P01569

Target Background

Produced by macrophages, IFN-alpha have antiviral activities. Interferon stimulates the production of two enzymes: a protein kinase and an oligoadenylate synthetase.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 22 kDa

Positive control: 231 cell, Human breast cancer tissue, Raji cell, Human fetal liver tissue lysates

Recommended dilution: 500-2000

Gel: 12% SDS-PAGE


Lysate: 40 µg

Lane 1-4: 231 cell, Human breast cancer tissue, Raji cell,
Human fetal liver tissue lysates

Primary antibody: ml164234 (IFNA5 Antibody) at dilution
1/400

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 90 seconds



ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn