

Anti-PONC antibody

Cat. No.	ml161037
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-PONC rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human PONC
Reactivity	Human, Mouse, Rat
Content	0.8 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	PONC
Full name	prepronociceptin
Synonyms	PPNOC
Swissprot	Q13519

Target Background

Nociceptin is the ligand of the opioid receptor-like receptor (OPRL1). It may act as a transmitter in the brain by modulating nociceptive and locomotor behavior. May be involved in neuronal differentiation and development

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 20 kDa

Positive control: Human liver cancer tissue

Recommended dilution: 1000-5000

Gel: 10% SDS-PAGE


Lysate: 30 µg

Lane: Human liver cancer tissue

Primary antibody: ml161037 (PONC Antibody) at dilution 1/1200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 second



ELISA

Recommended dilution: 2000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn