

Anti-KHSRP antibody

Cat. No.	ml164352
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-KHSRP rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human KHSRP
Reactivity	Human, Mouse, Rat
Content	1.38 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	KHSRP
Full name	KH-type splicing regulatory protein
Synonyms	FBP2; KSRP; FUBP2
Swissprot	Q92945

Target Background

The KHSRP gene encodes a multifunctional RNA-binding protein implicated in a variety of cellular processes, including transcription, alternative pre-mRNA splicing, and mRNA localization.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 73 kDa

Positive control: HepG2 and Hela cell lysates

Recommended dilution: 1000-5000

Gel: 8%SDS-PAGE


Lysate: 40 µg

Lane: HepG2 and Hela cell lysates

Primary antibody: ml164352(KHSRP Antibody) at dilution 1/1000

Secondary antibody: Goat anti rabbit IgG at 1/5000 dilution

Exposure time: 6 minutes



ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn