

Anti-GST Tag antibody

 Cat. No.
 ml190101

 Package
 25 μl/100 μl/200 μl

 Storage
 -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview Description Applications Immunogen Reactivity Content Host species Ig class Purification

Anti-GST Tag mouse monoclonal antibody ELISA, WB Full length fusion protein Fusion proteins containing the GST Tag 1.5 mg/ml Mouse Immunogen-specific mouse IgG Antigen affinity purification

Target information Symbol Full name Synonyms Swissprot

GST Tag Glutathione S-Transferase

Target Background


GST (Glutathione S-Transferase) is a 26kDa protein which was widely used in the pGEX family of GST plasmid expression vectors as a fusion protein with foreign proteins.GST isoenzymes appear to play a central role in the parasite detoxification system. Other functions are also suspected including a role in increasing the solubility of haematin in the parasite gut.



订购热线: 4008-898-798

Applications Western blotting Predicted band size:55 kDa Positive control:Recombinant GST tag protein Recommended dilution: 1000-5000

Gel: 10%SDS-PAGE Lysate: 2 µg Lane: Recombinant GST tag protein Primary antibody: ml190101(GST Tag Antibody) at dilution 1/1000 Secondary antibody: Goat anti mouse IgG at 1/8000 dilution Exposure time: 10 seconds



ELISA

Recommended dilution:1000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

- 邮箱: mlbio_cn@yeah.net
- 网址: www.mlbio.cn