

Anti-HSP90AA1 antibody

Cat. No.	ml220009
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-HSP90AA1 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human HSP90AA1
Reactivity	Human, Mouse, Rat
Content	0.6 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	HSP90AA1
Full name	heat shock protein 90kDa alpha (cytosolic), class A member 1
Synonyms	EL52; HSPN; LAP2; HSP86; HSPC1; HSPCA; Hsp89; Hsp90; LAP-2; HSP89A; HSP90A; HSP90N; HSPCAL1; HSPCAL4
Swissprot	P07900

Target Background

The protein encoded by this gene is an inducible molecular chaperone that functions as a homodimer. The encoded protein aids in the proper folding of specific target proteins by use of an ATPase activity that is modulated by co-chaperones. Two transcript variants encoding different isoforms have been found for this gene

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 85 kDa

Positive control: A549 and 293T cells, mouse brain tissue, K562 and Hela cells

Recommended dilution: 500-2000

Gel: 6%SDS-PAGE


Lysate: 40 µg

Lane 1-5: A549 cells, 293T cells, mouse brain tissue, K562 cells, Hela cells

Primary antibody: ml220009(HSP90AA1 Antibody) at dilution 1/350

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 5 seconds



ELISA

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn