

Anti-ACPP antibody

Cat. No.	ml220043
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-ACPP rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human ACPP
Reactivity	Human, Mouse
Content	0.3 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	ACPP
Full name	acid phosphatase, prostate
Synonyms	PAP, ACP3, 5'-NT, ACP-3
Swissprot	P15309

Target Background

This gene encodes an enzyme that catalyzes the conversion of orthophosphoric monoester to alcohol and orthophosphate. It is synthesized under androgen regulation and is secreted by the epithelial cells of the prostate gland. An alternatively spliced transcript variant encoding a longer isoform has been found for this gene. This isoform contains a transmembrane domain and is localized in the plasma membrane-endosomal-lysosomal pathway.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size:45 kDa

Positive control:Mouse liver tissue

Recommended dilution: 200-1000

Gel: 10%SDS-PAGE


Lysate: 40 μ g

Lane: Mouse liver tissue

Primary antibody: ml220043(ACPP Antibody) at dilution 1/250

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 10 seconds



ELISA

Recommended dilution:1000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn