

Anti-BAK1 antibody

Cat. No.	ml220146
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-BAK1 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human BAK1
Reactivity	Human
Content	0.3 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	BAK1
Full name	BCL2-antagonist/killer 1
Synonyms	BAK; CDN1; BCL2L7; BAK-LIKE
Swissprot	Q16611

Target Background

The protein encoded by this gene belongs to the BCL2 protein family. BCL2 family members form oligomers or heterodimers and act as anti- or pro-apoptotic regulators that are involved in a wide variety of cellular activities. This protein localizes to mitochondria, and functions to induce apoptosis. It interacts with and accelerates the opening of the mitochondrial voltage-dependent anion channel, which leads to a loss in membrane potential and the release of cytochrome c. This protein also interacts with the tumor suppressor P53 after exposure to cell stress.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 23 kDa

Positive control: Hela, A431, 293T and Jurkat cells

Recommended dilution: 500-2000

Gel: 12% SDS-PAGE


Lysate: 40 µg

Lane 1-4: Hela cells, A431 cells, 293T cells, Jurkat cells

Primary antibody: ml220146(BAK1 Antibody) at dilution 1/550

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 5 seconds



ELISA

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn