

Anti-APOBEC3C antibody

Cat. No.	ml161295
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-APOBEC3C rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human APOBEC3C
Reactivity	Human
Content	0.3 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	APOBEC3C
Full name	apolipoprotein B mRNA editing enzyme, catalytic polypeptide-like 3C
Synonyms	A3C; PBI; ARP5; ARDC2; ARDC4; APOBEC1L; bK150C2.3
Swissprot	Q9NRW3

Target Background

This gene is a member of the cytidine deaminase gene family. It is one of seven related genes or pseudogenes found in a cluster thought to result from gene duplication, on chromosome 22. Members of the cluster encode proteins that are structurally and functionally related to the C to U RNA-editing cytidine deaminase APOBEC1. It is thought that the proteins may be RNA editing enzymes and have roles in growth or cell cycle control.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 23 kDa

Positive control: PC3 cells

Recommended dilution: 200-1000

Gel: 12% SDS-PAGE

Lysate: 40 µg

Lane: PC3 cells

Primary antibody: ml161295 (APOBEC3C Antibody) at
dilution 1/350

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute



ELISA

Recommended dilution: 1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn