

Anti-FAM96A antibody

Cat. No.	ml223839
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-FAM96A rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Full length fusion protein
Reactivity	Human, Mouse
Content	0.6 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	FAM96A
Full name	family with sequence similarity 96, member A
Synonyms	
Swissprot	Q9H5X1

Target Background

May play a role in chromosome segregation through establishment of sister chromatid cohesion.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 18 kDa

Positive control: Human fetal intestines tissue

Recommended dilution: 200-1000

Gel: 12%SDS-PAGE


Lysate: 40 µg

Lane: Human fetal intestines tissue

Primary antibody: ml223839(FAM96A Antibody) at dilution 1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 30 seconds



ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn