

Anti-AHRR antibody

Cat. No.	ml161444
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-AHRR rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human AHRR
Reactivity	Human, Mouse
Content	1.5 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	AHRR
Full name	aryl-hydrocarbon receptor repressor
Synonyms	AHH; AHHR; bHLHe77
Swissprot	A9YTQ3

Target Background

The protein encoded by this gene participates in the aryl hydrocarbon receptor (AhR) signaling cascade, which mediates dioxin toxicity, and is involved in regulation of cell growth and differentiation. It functions as a feedback modulator by repressing AhR-dependent gene expression. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 76 kDa

Positive control: Mouse liver tissue and HepG2 cells

Recommended dilution: 500-2000

Gel: 6%SDS-PAGE


Lysate: 40 µg

Lane 1-2: Mouse liver tissue, HepG2 cells

Primary antibody: ml161444(AHRR Antibody) at dilution 1/430

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 2 minutes



ELISA

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn