

## Anti-CBX6 antibody

 Cat. No.
 ml220364

 Package
 25 μl/100 μl/200 μl

 Storage
 -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview	
Description	Anti-CBX6 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human CBX6
Reactivity	Human
Content	0.4 mg/ml
Host species	Rabbit
lg class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification
Target information	
Symbol	CBX6
Full name	chromobox homolog 6

O95503

## **Target Background**

Synonyms Swissprot


CBX6 (chromobox protein homolog 6) is a 412 amino acid protein that is a component of the chromatin-associated PcG complex. Through interaction with PcG, CBX6 modifies chromatin and maintains the transcriptionally repressive state of genes. Members of the chromobox domain protein family are characterized by the presence of a 37 amino acid chromobox (CHRromatin Organization MOdifier) domain. This domain plays a mechanistic role in targeting chromodomain proteins to specific regions of the nucleus.



订购热线: 4008-898-798

Applications Western blotting Predicted band size:44 kDa Positive control:Hela and 293T cells Recommended dilution: 500-2000

Gel: 10%SDS-PAGE Lysate: 40 µg Lane 1-2: Hela cells, 293T cells Primary antibody: ml220364(CBX6 Antibody) at dilution 1/600 Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution Exposure time: 40 seconds



## ELISA

Recommended dilution:2000-5000

Wit produce 联系电话: 4008-898-798, 021-61725725

## 联系QQ: 2881505695, 2881505696

- 邮箱: mlbio\_cn@yeah.net
- 网址: www.mlbio.cn