

# Anti-MCM8 antibody

Cat. No. ml220967 Package 25 µl/100 µl/200 µl -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol Storage

Product overview	
Description	Anti-MCM8 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human MCM8
Reactivity	Human, Mouse, Rat
Content	0.3 mg/ml
Host species	Rabbit
lg class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification
Target information	
Symbol	MCM8
Full name	minichromosome maintenance complex component 8
Synonyms	REC, C20orf154, dJ967N21.5
Swissprot	Q9UJA3
	certs
Target Background	

### **Target Background**

The protein encoded by this gene is one of the highly conserved mini-chromosome maintenance proteins (MCM) that are essential for the initiation of eukaryotic genome replication. The hexameric protein complex formed by the MCM proteins is a key component of the pre-replication complex (pre\_RC) and may be involved in the formation of replication forks and in the recruitment of other DNA replication related proteins. This protein contains the central domain that is conserved among the MCM proteins. This protein has been shown to co-immunoprecipitate with MCM4, 6 and 7, which suggests that it may interact with other MCM proteins and play a role in DNA replication. Alternatively spliced transcript variants encoding distinct isoforms have been described.



## 订购热线: 4008-898-798

Applications Western blotting Predicted band size:94 kDa Positive control:Jurkat cells Recommended dilution: 200-1000

Gel: 10%SDS-PAGE Lysate: 40 µg Lane: Jurkat cells Primary antibody: ml220967(MCM8 Antibody) at dilution 1/300 Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution Exposure time: 1 minute



### ELISA

Recommended dilution:1000-2000

Wit produce 联系电话: 4008-898-798, 021-61725725

## 联系QQ: 2881505695, 2881505696

- 邮箱: mlbio\_cn@yeah.net
- 网址: www.mlbio.cn