

Anti-DTX3L antibody

Cat. No.	ml220548
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-DTX3L rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human DTX3L
Reactivity	Human
Content	0.5 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	DTX3L
Full name	deltex 3 like, E3 ubiquitin ligase
Synonyms	BBAP
Swissprot	Q8TDB6

Target Background

The Deltex family is responsible for influencing Notch signaling and may regulate transcription through interactions with specific transcription factors. Deltex proteins have a basic N-terminus, a central proline-rich region and a C-terminal RING-type zinc finger domain, a motif often found in ubiquitin-protein isopeptide ligases (E3). The RING-type zinc finger domain binds two Zn²⁺ atoms and forms a cross-brace motif that is essential for many proteins involved the ubiquitination pathway. DTX3L (Deltex-3-like), also known as BBAP, is a 740 amino acid protein that is similar to Deltex-3 and acts as a ubiquitin ligase in vitro. DTX3L can heterodimerize with Deltex-1, a transcriptional regulator, thereby enhancing the activity of the E3 ubiquitin ligase complex and increasing the influence of E3 on the Notch signaling pathway.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 84 kDa

Positive control: A431 cell lysate

Recommended dilution: 200-1000

Gel: 6%SDS-PAGE


Lysate: 40 µg

Lane: A431 cell lysate

Primary antibody: ml220548(DTX3L Antibody) at dilution 1/300

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 5 seconds



ELISA

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn