

Anti-ADAM21 antibody

Cat. No.	ml221508
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-ADAM21 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human ADAM21
Reactivity	Human, Mouse
Content	0.36 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	ADAM21
Full name	ADAM metalloproteinase domain 21
Synonyms	ADAM31; ADAM 21
Swissprot	Q9UKJ8

Target Background

This gene encodes a member of the ADAM (a disintegrin and metalloprotease domain) family. Members of this family are membrane-anchored proteins structurally related to snake venom disintegrins, and have been implicated in a variety of biological processes involving cell-cell and cell-matrix interactions, including fertilization, muscle development, and neurogenesis. The expression of this gene is testis-specific.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size:81 kDa

Positive control:Mouse kidney tissue lysate

Recommended dilution: 500-2000

Gel: 8%SDS-PAGE


Lysate: 40 μ g

Lane: Mouse kidney tissue lysate

Primary antibody: ml221508(ADAM21 Antibody) at dilution 1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 40 seconds



ELISA

Recommended dilution:5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn