

订购热线: 4008-898-798

Anti-GAS2 antibody

Cat. No.	ml220674
Package	25 μl/100 μl/200 μl
Storage	-20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol
Product overview	
Description	Anti-GAS2 rabbit polyclonal antibody
Applications	ELISA, WB, IHC
Immunogen	Full length fusion protein
Reactivity	Human, Mouse
Content	0.3 mg/ml
Host species	Rabbit
lg class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification
Target information	prov
Symbol	GAS2
Full name	growth arrest-specific 2
C	00



Synonyms

Swissprot 043903

Target Background

The protein encoded by this gene is a caspase-3 substrate that plays a role in regulating microfilament and cell shape changes during apoptosis. It can also modulate cell susceptibility to p53-dependent apoptosis by inhibiting calpain activity. Multiple alternatively spliced variants, encoding the same protein, have been identified.




Applications

Immunohistochemistry

Predicted cell location: cytoplasm

Positive control: Human thyroid cancer

Recommended dilution: 25-100



The image on the left is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using ml220674(GAS2 Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: ×200)

Good elisa kit producert

Western blotting

Predicted band size:35 kDa

Positive control:Human fetal liver tissue

Recommended dilution: 200-1000



Gel: 12%SDS-PAGE	kDa
Lysate: 40 µg	95 — 72 — 55 —
Lane: Human fetal liver tissue	36 —
Primary antibody: ml220674(GAS2 Antibody) at dilution 1/200	28—
Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution	17 —
Exposure time: 5 seconds	10—

ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn