

订购热线: 4008-898-798

Anti-DIRAS3 antibody

Cat. No. ml221695

Package 25 μl/100 μl/200 μl

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-DIRAS3 rabbit polyclonal antibody

Applications ELISA, WB

Immunogen Fusion protein of human DIRAS3

ReactivityHuman, RatContent0.6 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol DIRAS3

Full name DIRAS family GTPase 3

Synonyms ARHI; NOEY2
Swissprot O95661

Target Background

This gene encodes a member of the ras superfamily. This gene is imprinted gene with monoallelic expression of the paternal allele which is associated with growth suppression. The encoded protein acts as a tumor suppressor whose function is abrogated in many ovarian and breast cancers. This protein may also play a role autophagy in certain cancer cells by regulating the autophagosome initiation complex.



订购热线: 4008-898-798

Applications

Western blotting

Predicted band size:26 kDa

Positive control:Rat brain tissue and Rat liver tissue lysates

Recommended dilution: 200-1000

Gel: 12%SDS-PAGE

Lysate: 40 µg


Lane 1-2: Rat brain tissue and Rat liver tissue lysates

Primary antibody: ml221695(DIRAS3 Antibody) at dilution

1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 3 minutes



ELISA

Recommended dilution:5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn