

订购热线: 4008-898-798

Anti-EZH2 antibody

 Cat. No.
 ml221769

 Package
 25 μl/100 μl/200 μl

 Storage
 -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Dreduct even deve	
Product overview	
Description	Anti-EZH2 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human EZH2
Reactivity	Human, Mouse
Content	0.2 mg/ml
Host species	Rabbit
lg class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification
Target information	
Symbol	EZH2
Full name	enhancer of zeste homolog 2 (Drosophila)
Synonyms	WVS; ENX1; EZH1; KMT6; WVS2; ENX-1; EZH2b; KMT6A

Q15910

Target Background

Swissprot


This gene encodes a member of the Polycomb-group (PcG) family. PcG family members form multimeric protein complexes, which are involved in maintaining the transcriptional repressive state of genes over successive cell generations. This protein associates with the embryonic ectoderm development protein, the VAV1 oncoprotein, and the X-linked nuclear protein. This protein may play a role in the hematopoietic and central nervous systems. Multiple alternatively splcied transcript variants encoding distinct isoforms have been identified for this gene.



订购热线: 4008-898-798

Applications Western blotting Predicted band size:85 kDa Positive control:K562 cells Recommended dilution: 500-2000

Gel: 8%SDS-PAGE Lysate: 40 µg Lane: K562 cells Primary antibody: ml221769(EZH2 Antibody) at dilution 1/300 Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution Exposure time: 15 seconds



ELISA

Recommended dilution:2000-5000

Wit produce 联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

- 邮箱: mlbio_cn@yeah.net
- 网址: www.mlbio.cn