

## Anti-NOL9 antibody

<b>Cat. No.</b>	ml221873
<b>Package</b>	25 µl/100 µl/200 µl
<b>Storage</b>	-20°C, pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol

### Product overview

<b>Description</b>	Anti-NOL9 rabbit polyclonal antibody
<b>Applications</b>	ELISA, WB
<b>Immunogen</b>	Fusion protein of human NOL9
<b>Reactivity</b>	Human
<b>Content</b>	0.3 mg/ml
<b>Host species</b>	Rabbit
<b>Ig class</b>	Immunogen-specific rabbit IgG
<b>Purification</b>	Antigen affinity purification

### Target information

<b>Symbol</b>	NOL9
<b>Full name</b>	nucleolar protein 9
<b>Synonyms</b>	Grc3; NET6
<b>Swissprot</b>	Q5SY16

### Target Background

Nucleolar protein 9 is a protein that in humans is encoded by the NOL9 gene. Polynucleotide 5'-kinase involved in rRNA processing. The kinase activity is required for the processing of the 32S precursor into 5.8S and 28S rRNAs, more specifically for the generation of the major 5.8S(S) form. In vitro, has both DNA and RNA 5'-kinase activities. Probably binds RNA.

订购热线: 4008-898-798

**Applications**

**Western blotting**

Predicted band size:79 kDa

Positive control:Hela cells

Recommended dilution: 500-2000

Gel: 10%SDS-PAGE

Lysate: 40 µg

Lane: Hela cells

Primary antibody: ml221873(NOL9 Antibody) at dilution 1/1250

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 5 seconds



**ELISA**

Recommended dilution:2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio\_cn@yeah.net

网址: www.mlbio.cn