

## Anti-KHK antibody

 Cat. No.
 ml221887

 Package
 25 μl/100 μl/200 μl

 Storage
 -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview	
Description	Anti-KHK rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human KHK
Reactivity	Human, Mouse
Content	0.4 mg/ml
Host species	Rabbit
lg class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification
Target information	
Symbol	КНК
Full name	ketohexokinase (fructokinase)
Synonyms	

P50053

## **Target Background**

Swissprot


This gene encodes ketohexokinase that catalyzes conversion of fructose to fructose-1-phosphate. The product of this gene is the first enzyme with a specialized pathway that catabolizes dietary fructose. Alternatively spliced transcript variants encoding different isoforms have been identified.



订购热线: 4008-898-798

Applications Western blotting Predicted band size:33 kDa Positive control:Mouse liver and kidney tissue, human fetal liver tissue Recommended dilution: 1000-5000

Gel: 15%SDS-PAGE Lysate: 40 µg Lane 1-3: Mouse liver tissue, Mouse kidney tissue, human fetal liver tissue Primary antibody: ml221887(KHK Antibody) at dilution 1/1350 Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution Exposure time: 3 seconds



## ELISA

Recommended dilution:2000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio\_cn@yeah.net

网址: www.mlbio.cn