

Anti-BRD9 antibody

Cat. No.	ml224272
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-BRD9 rabbit polyclonal antibody
Applications	ELISA, WB, IHC
Immunogen	Fusion protein of human BRD9
Reactivity	Human, Mouse
Content	0.8 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	BRD9
Full name	bromodomain containing 9

Synonyms PRO9856; LAVS3040

Swissprot Q9H8M2

Target Background

May play a role in chromatin remodeling and regulation of transcription.


Applications

Immunohistochemistry

Predicted cell location: Cell membrane

Positive control: Human liver cancer

Recommended dilution: 50-300



The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using ml224272(BRD9 Antibody) at dilution 1/80, on the right is treated with fusion protein.
(Original magnification: ×200)

Western blotting

Predicted band size: 67 kDa

Positive control: HEPG2, HeLa, 293T cell lysates

Recommended dilution: 500-2000

Gel: 8%SDS-PAGE


Lysate: 40 μ g

Lane 1-3: HEPG2, Hela, 293T cell lysates

Primary antibody: ml224272(BRD9 Antibody) at dilution 1/550

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 5 seconds



ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn