



订购热线: 4008-898-798

## Anti-BCAM antibody

Cat. No.	ml222068
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

### Product overview

Description	Anti-BCAM rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human BCAM
Reactivity	Human
Content	0.3 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

### Target information

Symbol	BCAM
Full name	basal cell adhesion molecule (Lutheran blood group)
Synonyms	AU; LU; CD239; MSK19
Swissprot	P50895

### Target Background

This gene encodes Lutheran blood group glycoprotein, a member of the immunoglobulin superfamily and a receptor for the extracellular matrix protein, laminin. The protein contains five extracellular immunoglobulin domains, a single transmembrane domain, and a short C-terminal cytoplasmic tail. This protein may play a role in epithelial cell cancer and in vaso-occlusion of red blood cells in sickle cell disease. Polymorphisms in this gene define some of the antigens in the Lutheran system and also the Auberger system. Inactivating variants of this gene result in the recessive Lutheran null phenotype, Lu(a-b-), of the Lutheran blood group. Two transcript variants encoding different isoforms have been found for this gene.



Good elisakit producers

订购热线: 4008-898-798

#### Applications

##### Western blotting

Predicted band size: 67 kDa

Positive control: 293T and 231 cells

Recommended dilution: 200-1000

Gel: 6% SDS-PAGE


Lysate: 40 µg

Lane 1-2: 293T cells, 231 cells

Primary antibody: ml222068(BCAM Antibody) at dilution 1/400

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute



#### ELISA

Recommended dilution: 1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio\_cn@yeah.net

网址: www.mlbio.cn